

The image shows a white, handheld URA Reader device with a circular sensor area on the front. The device has a sleek, modern design with rounded edges and a white finish. The brand name 'QASSAY' is visible on the top part of the device. The background is a light teal color with abstract, circular patterns and a large, faint 'Q' shape. The overall aesthetic is clean and professional, suggesting a high-tech medical or laboratory device.

URA Reader

Welcome to a new era in Urinalysis

APPLICATION	Urine Strip Analyzer
INTERFACE	Qassay Smartphone App (Android / iOS)
DIMENSIONS	170 mm * 50 mm * 16 mm
WEIGHT	59 g
NUMBER OF PADS	12-13
OPERATING MODE	Cassette manual extraction
SENSOR	High resolution 10 channel spectrometer (350 - 1050 nm).
MEASUREMENT TIME	Approx. 2s.
POWER SUPPLY	USB-C Rechargeable Li-Ion battery.
DATA MANAGEMENT	Web panel for data management (API Rest)
CONNECTIVITY	Bluetooth 5.0
TEST IDENTIFICATION	QR code read through smartphone application.
OPTIMAL CONDITIONS	Between 5°C and 40°C; and <90% RH



QASSAY

QASSAY.COM